

GAMMEX[®] Non-Latex PI Ortho

S u r g i c a l G l o v e s

Polyisoprene, skin-friendly, ultra-thick surgical glove to withstand the rigors of surgery, especially where power tools and bone fragments are present

A next-generation, skin-friendly polyisoprene (PI) glove free of chemical accelerators known to cause Type IV glove allergies and sensitivities.* Offering 50% extra thickness on average, this non-latex glove will help withstand the rigors of high risk, trauma and reconstructive procedures while offering exceptional grip and comfort.

KEY FEATURES AND BENEFITS:

- 50% thicker on average than Ansell's standard powder-free latex surgical gloves for outstanding protection
- Free from chemical accelerators known to cause chemical Type IV allergies and sensitivities for optimized allergenic profile
- Unique, anti-slip external surface treatment for exceptional handling and grip
- The DERMASHIELD™ Technology internal coating facilitates dry and damp donning as well as intra-operative changes

**Data on file: HRIPT (Human Repeat Insult Patch Testing) conducted by an independent laboratory in the U.S. PI products featuring PI-KARE Technology may not yet be available in your country. Please contact your Ansell representative.*

Recommended For

- All surgeries, especially high risk procedures in orthopaedics, trauma, and cardiac
- Damp donning
- Protection from chemical accelerators such as thiurams, carbamates and mercaptos known to cause chemical Type IV allergies

FEATURED TECHNOLOGIES



GAMMEX[®] Innovation that Protects

Driven by innovation, GAMMEX[®] surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

GAMMEX[®] Non-Latex PI Ortho

Surgical Gloves

PRODUCT INFORMATION

Material	PI (Polyisoprene)
Powder Content	Powder-Free
Color	Green
Thickness	Ortho
Cuff Style	Beaded with SUREFIT™ Technology
External Glove Surface	Smooth with Micro-Textured Finish
Internal Glove Surface	Polymer Coated with DERMASHIELD™ Technology
Allergy Prevention	Latex (Type I)
Grip Level	Moderate
Double Gloving Recommendation	No
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)

ORDERING INFORMATION

Size	6	6.5	7	7.5	8	8.5	9
Ansell Product Code	20686560	20686565	20686570	20686575	20686580	20686585	20686590

PACKAGING AND STORAGE

Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition

STANDARDS AND CERTIFICATIONS

Product Standards	AS/NZS 4179, ASTM D3577, ASTM D6978, ASTM D7160, EN 420:2003 + A1:2009
Quality / Environmental Standards	EN556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email: insidesalesus@ansell.com

Central & South America
Tel: +52(442) 296 20 50
Email: csllac@ansell.com
Brazil
Tel: +55-11-3356-3100
Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com