

#whdentalwerk



video.wh.com

Now at
your dealer
or **wh.com**



The formula for a long instrument life



f1 service
oil

The More
Intelligent Formula

Prolong the lifespan of your instruments

For Every Application

W&H Service Oil F1 is a product of W&H research. It has been developed hand in hand with modern W&H instrument technology. This oil gives you the confidence of knowing that you are providing your W&H instruments with the best possible care:

- › Turbines
- › Straight and contra-angle handpieces
- › Surgical instruments
- › Air scalers
- › Air motors



W&H dental drives achieve top speeds of up to 390,000 r.p.m. F1 provides excellent lubrication and internal cleaning thereby ensuring fault-free operation of the instruments.



Perfect for manual maintenance
W&H Service Oil F1 in the MD-400 spray can.



Perfect in the Assistina
W&H Service Oil F1 in the MD-500 bottle for use in the Assistina.

With the test champion oil for all standard drives

The More Intelligent Service Oil
W&H Service Oil F1 is an special oil of maximum purity, which has been specially developed for the demanding and extreme requirements of dental instruments. Compared with other well-known brand oils.


Sterilization resistant 135°C
155
F1 is of course fully sterilizable, absolutely resistant to hydrolysis, stable against oxidation and is shown to be highly compatible with disinfectants.

Physiologically o.k.
F1 contains no toxic constituents, is physiologically harmless and environmentally compatible. The substances it contains comply with the applicable international regulations.

High material compatibility
F1 also safeguards the durability of plastic parts, coatings, seals, glass etc.

Enormously economical
Your instruments will thank you for regular maintenance with F1 by requiring fewer repairs. The advantages: considerable savings in cost and increased profitability of your W&H instruments.

Short lubrication intervals – long service life
Intensive maintenance will considerably increase the life of your instruments.



Intelligent oil + short maintenance intervals = longer lifespan

Oil spray can
Type MD-400
Contents 400 ml

Cartridge for Assistina
3x3 / 3x2
Type MD-400
Contents 200 ml

Bottle for Assistina
Type MD-500
Contents 500 ml



Order your Service Oil F1 through your dental supplier.
Suitable spray caps can be found on the Internet at **wh.com**

W&H Dentalwerk Bürmoos GmbH
Ignaz-Glaser-Straße 53, Postfach 1
5111 Bürmoos, **Austria**
t +43 6274 6236-0
f +43 6274 6236-55
office@wh.com
wh.com

Pictures are for illustrative purposes only. Additional equipment and accessories shown are not included as standard.

20337 AEN Rev. 005 / 26.07.2019
Subject to alterations