


GRIP-LOK[®]



Catheter and Line Securement
Safe, Simple, Secure

GRIP-LOK®

Grip-Lok® is a stabilization and securement device that is strong enough to lock catheters, tubes and lines securely in place, yet versatile enough for almost any securement situation. Engineered Stabilization Devices (ESD) are recommended by INS to promote consistent practice, reduce vascular access device (VAD) motion and complications, and reduce interruption of needed infusion therapy.**

** "Infusion Therapy Standards of Practice." Journal of Infusion Nursing 39:IS (2016): S72-74. Web.

Safe...

Grip-Lok® is packaged sterile, allowing you to use Grip-Lok® *close to the insertion site or right in the sterile field.*

Simple...

Grip-Lok® is applied in the same simple manner as a bandage and can be easily placed and adjusted.

Secure...

Durably engineered, Grip-Lok® helps prevent dislodgement and migration of the catheter. This can help reduce lost lines and patient discomfort.

Contents

General Securement	3
Peripheral, Arterial & Fistula	4
CVC & PICC.	5
Foley	6
Nasal Gastric	7



General Securement



Small

Grip-Lok® Small securement device is strong enough to lock lines and tubes 1.5-4.5mm in diameter.

Overall length: 3" (7.6cm)

Medium

Grip-Lok® Medium securement device is versatile, strong, and can hold lines and tubes 3-8mm in diameter.

Overall length: 3.5" (8.9cm)



Large

Grip-Lok® Large securement device is strong and contains a large footprint, and can hold lines and tubes 6-13mm in diameter.

Overall length: 6.5" (16.5cm)

PRODUCT #	DESCRIPTION	U/M
3200S	Small Securement Device	100/bx
3300M	Medium Securement Device	100/bx
3400L	Large Securement Device	100/bx

Peripheral & Arterial Catheter Securement

Peripheral IV and Arterial Catheter Securement

Grip-Lok® Peripheral and Arterial securement devices secure the catheter and Luer hub with a unique wrap around design. The device provides securement and protection for the skin from hard plastic parts.

Overall length: 3.5" (8.9cm)



Dialysis Fistula Needle

Designed to secure Dialysis Fistula Needles.

Overall length: 4.5" (11.4cm)
Two devices per sterile pouch.
One box contains 100 sterile pouches.

PRODUCT #	DESCRIPTION	U/M
3300MIV	Peripheral IV and Arterial Catheter Securement Device	100/bx
3300MART	Arterial and Peripheral IV Catheter Securement Device	100/bx
3300MIV-KIT-2	IV Start Kit with 1-Grip-Lok® 3300MIV, 1 Tape Strip, 1-2.37" x 2.75" Transparent Dressing	100/bx
3500DFN	Dialysis Fistula Needle Securement Device	100/bx*

* Two securement devices per sterile package.



3300MIV and 3300MART will secure ported catheters.



3300MIV and 3300MART will secure needleless connectors.

CVC & PICC Securement

CVC & PICC Securement



Medium Wide Adhesive

Grip-Lok® Medium Wide Adhesive will adhere and secure catheters made from silicone, PVC and other plastic compounds.

Overall length: 3.5" (8.9cm)

PICC and CVC Securement

Grip-Lok® 3601CVC for PICC and CVC securement provides an adhesive area that will stick to a wide range of hubs and delivers protection from accidental line pulls.

Overall length: 3.5" (8.9cm)

The 3601CVC can secure a catheter on the chest or neck and is designed to provide stable suture-free securement at the catheter attachment points.



Extra Small PICC and Line Securement

Grip-Lok® Extra Small PICC securement device will fit a range of hubs made from silicone and other plastic compounds. Additionally, this device contains zinc oxide skin contact adhesive.

Designed for neonatal sized PICC lines. | Overall length: 2.25" (5.7cm)

PRODUCT #	DESCRIPTION	U/M
3300MWA	Medium Wide Securement Device	100/bx
3601CVC	PICC and CVC Securement Device	100/bx
2200NUZA	Extra Small PICC and Line Securement Device	100/bx

Foley Catheter Securement



Grip-Lok® Large Foley Catheter Securement

Grip-Lok® for Foley Catheter does not utilize adhesive to secure the catheter in place, allowing for freedom of movement and patient comfort. Grip-Lok® Large Foley Catheter Securement is designed to hold the catheter in a channel created by a foam insert and a fabric loop over-strap. The degree of rotational and horizontal movement can be adjusted using the top fabric over-straps.

Designed for tubes sized 12-30 French in diameter.
Overall length: 6.25" (15.9cm)

Small Foley Catheter Securement

Grip-Lok® Small Foley Catheter Securement works just like the larger 3400LFC, but the reduced size may be preferred for small patients.

Designed for tubes sized 12-30 French in diameter.
Overall length: 3.5" (8.9cm)



PRODUCT #	DESCRIPTION	U/M
3400LFC	Large Foley Securement Device	100/bx
3600PFC	Small Foley Catheter Securement Device	100/bx

Nasal Gastric Securement



Nasal Gastric Securement

Grip-Lok® Nasal Gastric offers firm securement of feeding tubes and will fit 3-5mm tube sizes.

Designed for tubes sized 1/8"-13/64" (3.0-5.0mm) in diameter.
Overall length: 1.75" (4.45cm)

Nasal Gastric Hydrocolloid

Grip-Lok® Nasal Gastric Hydrocolloid provides firm securement of feeding tubes and will fit tube sizes ranging from 1-2mm in diameter. With a very small footprint, the Nasal Gastric Hydrocolloid is designed to provide securement on the smallest of patients.

Designed for tubes sized 1/25"-5/64" (1.0-2.0mm) in diameter.
Overall length: 1.65" (4.2cm)



Large Nasal Gastric Securement

Grip-Lok® Large Nasal Gastric provides a wide adhesive area to distribute pull force across the skin surface.


Designed for tubes sized 1/16"-1/2" (1.5-13.0mm) in diameter.
Overall length: 4.25" (10.8cm)

PRODUCT #	DESCRIPTION	U/M
2100ANG	Nasal Gastric Securement Device	100/bx
2100NGH	Small Nasal Gastric Hydrocolloid Securement Device	50/bx
3401LNG	Large Nasal Gastric Securement Device	100/bx

Grip-Lok® products are available through your preferred distributor.

Product availability may vary by country/region. Please contact TIDI® Products for product, pricing and availability information in your area.

Visit www.tidiproducs.com/grip-lok for trial and inservice support and to download Product Brochures, Application Sheets, Instructions For Use, and Instructional Videos.



Distributed by:



Grip-Lok is a registered trademark of TIDI Products, LLC

HYPOALLERGENIC



Not made with natural rubber latex

EC REP Authorized European Representative

Blue Box Medical Limited
Unit 29, New Forest Enterprises Centre,
Rushington Business Park, Chapel Lane
Totton, Southampton, SO40 9LA, U.K.
Tel: 023 8066 9000; Fax: 023 8066 9002

